

UV LIGHT

Burton's Ultraviolet Light is the right light for a wide variety of applications.

Designed for numerous applications, Burton UV Lights are offered in models with magnifiers and white lights to best suit your needs. The UV light is also used in applications for criminology and industrial inspection.



PERFORMANCE RELIABILITY VALUE



Burton's UV lights are used worldwide in numerous applications.

The Right Light


MAGNIFIERS

UV LIGHT

The ideal light for numerous applications.

UV LIGHT FEATURES:

- > Special fluorescent bulbs produce UVA at approximately 360 nm wavelength
- > Available with magnifier and white bulbs
- > 4-light magnifier models include optically ground and polished glass 5 diopter magnifier lens with 2.25X total power magnification and focal length of 8"
- > 6,000 hour average bulb life
- > 115V model standard; 230V models available
- > 5-year warranty
- > Made in the USA



MODELS		
115V	230V	
#UV501	n/a	Hand-held UV Light with 2 UV Bulbs (115V)
#UV502	#UV522	UV Light with Magnifier and 2 UV Bulbs
#UV503	#UV523	UV Light with Magnifier and 2 White and 2 UV Bulbs

ACCESSORIES	
#1007090	Wall Bracket for UV Lights
#0001127PK	Set of 4 White Replacement Bulbs
#1003073PK	Set of 4 UV Replacement Bulbs



Luxo Medical Products

BURTON MEDICAL
 An ISO Certified Company
 21100 Lassen Street
 Chatsworth, CA 91311
 Phone: 800-444-9909, 818-701-8700
 Fax: 800-765-1770, 818-701-8725
www.burtonmedical.com