

burton

Luxo Medical Products

# XENALUX™

*Advanced design assures top performance and a noticeably brighter light beam.*

Powered by a 300-watt Xenon source and with light delivery via multi-fiber cables, Burton's XenaLux™ O.R. headlamp system offers superior intensity as well as remarkable reliability and value. The headlamp's optical system has a full range of articulation with a variable spot range, removable joystick adjustment and a new and improved vertically adjustable linkage arm for added reach and greater maneuverability. Unique design and optics allow the XenaLux™ Illuminator to provide superior features.



PERFORMANCE RELIABILITY VALUE



Burton's XenaLux™ 300-Watt Xenon Headlamp – powerful, reliable and lightweight!

*The Right Light*

MAJOR O.R.

# XENALUX™

The perfect complement to Burton's line of O.R. lights.

## XENALUX™ HEADLAMP FEATURES:

- > Extremely bright light of 120,000+ lux
- > Variable spot with flawless circular beam
- > Very light weight right angle headlamp
- > Preferred for comfort
- > Compact optical system
- > Fully-adjustable, uniquely-designed linkage arm allows for vertical adjustment of the optical module
- > Adapts to any Wolf port
- > Attractive standard 3.5 bundle cables – other cables including glass bifurcated, are available
- > Extraordinary warranty: 3 years

Virtually weightless headlamp minimizes neck and shoulder strain during lengthy surgeries – a noticeable difference in comfort.



Headband is fully adjustable and equipped with removable and replaceable sweat cushions.

The optical system has a full range of articulation with a variable spot range, removable baton adjustment and a new and improved vertically adjustable linkage arm for added reach and greater maneuverability.

Available in a unique breath-o-preme™, non-slip sweat-absorbent material for added comfort and convenience during long procedures.

## XENALUX™ ILLUMINATOR FEATURES:

- > Produces highest intensity white light of any 300-Watt Xenon source
- > Long bulb life of 500 hours
- > Quiet, cool venting system
- > Rotary Turret for Wolf, Storz, Olympus, ACMI and Pilling connectors
- > Extraordinary warranty: 3 years
- > Hour-in-use meter
- > Optional mobile floorstand offers vertical adjustment on upright assembly; the Illuminator's shelf includes a handle for easy positioning



Optional mobile stand for added convenience.



## MODELS

#XN10C 300-watt Illuminator with Headlamp

## COMPONENTS

#XN10H Headlamp System

#XN101 300-watt Illuminator

## ACCESSORIES

Call Burton for more information: 1-800-444-9909

**burton**

Luxo Medical Products

BURTON MEDICAL  
21100 Lassen Street  
Chatsworth, CA 91311  
Phone: 800-444-9909, 818-701-8700  
Fax: 800-765-1770, 818-701-8725  
www.burtonmedical.com