

LIGHT ACCESSORIES

Available for Outpatient® II, CoolSpot™ II, 12" (30 cm) Reflectors, Visionary®, Genie™ II, and Genesis® Plus lights. (Consult Customer Service on the appropriate accessories for each light.)



Burton Bulb Packs (see Bulb Matrix)

To ensure correct operation, prevent product failure and possible loss of warranty, use Burton halogen replacement bulbs. Burton bulbs are tested at our facilities to meet MIL-SPEC standards and are packed in custom boxes to signify they meet our performance requirements.



Disposable Light Handle Covers

Save time by eliminating the need to autoclave handles and allows compliance with OSHA standards for sterile procedures.



CoolSpot™ II Handle Mount Assembly

Select the CoolSpot™ II autoclavable handle for added maneuverability and convenience.