

# REFLECTORS

*Durable and powerful, Burton reflectors are ideal for both examinations and minor procedures.*

The “standard” in the medical industry, Burton Reflectors offer uncompromising quality and reliability. Boasting the best features and durability of any reflector-type light in its class, Burton Reflectors feature 100-watt halogen bulbs protected by UV/IR absorbent glass cylinders providing a cool light beam for added comfort during long hours of use. The Reflector is the ideal light for nearly all procedures.



PERFORMANCE RELIABILITY VALUE



Powerful, reliable illumination make Burton Reflectors the industry standard.

*The Right Light*



DIAGNOSTIC

# REFLECTORS

Powerful, reliable illumination make Burton Reflectors the industry standard.

## REFLECTOR FEATURES:

- > 3500 fc (38,000 Lux) at 24" (61 cm)
- > 3900° K Color Temperature
- > Includes a 100-watt halogen bulb with 2000-hour life
- > Bulb housed in a heat-absorbing, color-correcting glass cylinder
- > Proven design reduces shadows in field of vision
- > Parallel motion spring-tension mountings provide years of trouble-free use in the field
- > Battery option
- > Two SteriHandles and Adapters (available free with new units; can be purchased for older Reflectors)
- > 115V standard; 230V models also available
- > UL 60601-1/60601-2-41 Certified
- > 5-year warranty
- > Made in the USA



MODELS		
115V	230V	
#R400FL	#R402FL	12" Diameter Floorstand
#R400SC	#R402SC	12" Diameter Single Head Ceiling Mount

ACCESSORIES	
#1006760	Replacement Heat Glass Cylinder Kit
#1003016	29" Drop Tube Rod
#0006126PK	Set of 4 Halogen Reflector 100W Replacement Bulbs – 115V

NOTE: The optical data are nominal values based on measurements done according to the UL/IEC performance standards at 120V.



Luxo Medical Products

**BURTON MEDICAL**  
 An ISO Certified Company  
 21100 Lassen Street  
 Chatsworth, CA 91311  
 Phone: 800-444-9909, 818-701-8700  
 Fax: 800-765-1770, 818-701-8725  
[www.burtonmedical.com](http://www.burtonmedical.com)