

GLEAMER®

Aesthetics, advanced optics, and ergonomics assure uncompromising performance and comfort.

If your facility is searching for an affordable, all-purpose, rugged, compact halogen light, Burton's Gleamer® is your solution. The Gleamer® was designed to meet healthcare facilities' most demanding specifications for high intensity, cool operation, compact design, and value. Offering 4166 footcandles, a sturdy gooseneck, plus optional 16" (41cm) extension arm, and superior quality components, the Gleamer® provides high intensity, maximum flexibility and extraordinary reliability. Its unique venting system, energy-absorbing infrared exchanger and patented reflector-filter coatings provide a heat-free beam in a cool, yet compact housing.



PERFORMANCE RELIABILITY VALUE



Designed to meet even the most demanding healthcare facility's specifications for high intensity, cool operation, and low cost.

The Right Light

TASK & EXAMINATION

GLEAMER®

Burton's affordable, all-purpose halogen light.

GLEAMER® FEATURES:

- > 4166 footcandles (45,000 Lux) – Nearly 35% brighter!
- > 3700°K color temperature (natural coloring)
- > 35-watt halogen bulb; 4000-hour life
- > Easy bulb replacement without tools
- > Cool housing temperature and cool light beam; light can be in continuous use for hours without causing discomfort to physician or patient
- > Convenient positioning with ergonomic positioning ridges
- > Floorstand includes unique space-saving base (casters included)
- > High-quality, space-saving 4-caster base – safe and flexible portable movement with locking feature in 2 of 4 casters
- > 115V model standard; 230V model available
- > 5-year warranty
- > UL 60601-1/60601-2-41 Certified
- > Made in the USA




MODELS

STANDARD SPOT

115V	230V	
#GL30FL	#GL32FL	Gleamer® Floorstand (with casters)
#GL30W	#GL32W	Gleamer® with Extension Arm and Universal Wall Mount

WIDE BEAM SPOT

115V	230V	
#GLW30FL	#GLW32FL	Gleamer® Floorstand (with casters)
#GLW30W	#GLW32W	Gleamer® with Extension Arm and Universal Wall Mount



ACCESSORIES

#0009600PK Set of (4) 35-Watt Replacement Bulbs (for standard spot)

#0009608PK Set of (4) 35-Watt Replacement Bulbs (for wide beam spot)

burton

Luxo Medical Products

BURTON MEDICAL

An ISO Certified Company

21100 Lassen Street

Chatsworth, CA 91311

Phone: 800-444-9909, 818-701-8700

Fax: 800-765-1770, 818-701-8725

www.burtonmedical.com

NOTE: The optical data are nominal values based on measurements done according to the UL/IEC performance standards at 120V at 18".