

# LED 125

*The latest advancement in portable exam “HeadLite” technology.*

Burton’s LED 125 HeadLite is a durable, affordable alternative to similar filament bulb (Xenon, Halogen etc.) headlamps. Shock and vibration resistant, LEDs have no glass or filaments to break – making them much more durable than other light sources. The combination of a long LED life (approximately 50,000 hours of use) and small power needs (the LED 125 HeadLite will run continuously on a single set of batteries for over 100 hours) translates into lower operation, maintenance and upkeep costs for your practice. Also, the LED 125 HeadLite’s teal-colored light and precise beam is easier on the eyes than other headlamps.



PERFORMANCE RELIABILITY VALUE



The durable, affordable alternative to Xenon and Halogen headlamps.

*The Right Light*


TASK & EXAMINATION

# LED 125

Burton's affordable, reliable 1.25 watt LED exam headlight is almost maintenance-free.

## LED 125 HEADLITE FEATURES:

- > Crystal reflector tubes for extra brightness; approximately 743 footcandles (8,000 Lux)
- > Powder-coated metal casing protection
- > Nearly unbreakable Light Emitting Diode (LED)
- > Estimated 50,000-hour LED lifespan
- > Water resistant / shock proof
- > Weighs only 4.7 ounces
- > Includes a set of three 'AAA' batteries providing 100+ hours of life
- > Prism reflector system for optimal focus
- > 1-year warranty



85° hinged adjustable beam includes 10 positions for precise illumination

## MODELS

LED 125    LED 125 HeadLite 1.25 Watt

Secure, soft, no sweat, adjustable headband

Three AAA batteries deliver 100+ hours of operation



1.25 watt LED



**burton**

Luxo Medical Products

BURTON MEDICAL  
21100 Lassen Street  
Chatsworth, CA 91311  
USA

Phone: 800-444-9909, 818-701-8700  
Fax: 800-765-1770, 818-701-8725  
www.burtonmedical.com

NOTE: Product specifications may be modified by manufacturer at any time.