



Magnifiers

## Epic<sup>®</sup> Magnifier

The perfect combination of world-class optics, precision movement and a large viewing area.

Philips Burton's Epic<sup>®</sup> Magnifier commands attention based on aesthetics alone – but its real value lies in its world class optics, performance and precision movement. With its patented arm system, the Epic<sup>®</sup> is more rugged than any other magnifier. The arm stops exactly where positioned and does not require constant tightening or lock mechanism adjustments.

*The Right Light*

**PHILIPS**  
**burton**

## Features

- Parallel motion, patented enclosed arm system
- Quick start electronic ballast (flicker-free)
- Integral handles and cover for lens protection
- Available in floorstand and table/wall models



## Summary of technical data

Diopter of magnifying lens:	3.5
Viewing area diameter:	6.5" (16.5 cm)
Rated life of bulb:	10,000-hour average lamp life
Type of bulb:	22W
Certifications/approvals:	UL 153 US/CAN; CE marked
Product warranty:	1-year limited warranty



Note: The optical data are nominal values based on measurements calculated according to UL/IEC performance standards at 120V.

## Ordering Information

120V Model	230/240V Model	Description
EM60FL	EM62FL	Floorstand on Mobile Base (casters included)
EM60W	EM62W	Wall / Table Mount
Accessories		
#0009800PK		Epic® Magnifier Replacement Bulbs, 120V, 2/Box
#0009400PK		Epic® Magnifier Replacement Bulbs, 230/240V, 2/Box
#8970-001		4 Diopter Lens Attachment
#8970-002		6 Diopter Lens Attachment



Philips Burton  
21100 Lassen Street Chatsworth, CA 91311  
Phone: 818-701-8700 Fax: 818-701-8725  
www.burtonmedical.com

© 2011 Philips Burton. All rights reserved. Product designs protected by copyright.  
We reserve the right to change details of design, materials and finishes.

LIT-0007, Rev. 3