



# Veterinary Lighting

State-of-the-art lighting solutions for nearly all veterinary applications.

*The Right Light*

**PHILIPS**  
**burton**

# Veterinary Lighting

Uncompromising performance, style and value are guaranteed.

Philips Burton has just the right Veterinary lights to meet your needs, including energy and cost-efficient LED models. Offering Veterinary lights for nearly all applications, Philips Burton provides a range of products tailored to small animals as well as products specific for large animals. Our bright, cool lights lead in value and best of all, keep the veterinarian, staff and patient comfortable and safe.



Philips Burton provides a range of lighting products tailored for use with both small and large animals.

## The Right Light for Veterinary Applications

Application	Small Animal	Large Animal
General Task	LE-35	LE-35
Examination	Gleamer®, Super Exam 50® LED	Gleamer®, Super Exam 50® LED
Procedure	Outpatient® II, CoolSpot II®, AIM-50®	Outpatient® II, CoolSpot II®, AIM-50®
Minor Surgery	AIM-50®, AIM LED, CoolSpot II®	AIM-100®, AIM LED, CoolSpot II®
Major Surgery	AIM-100®, AIM-200® OR	AIM-200® OR, APEX
Magnifiers	Woods Magnifier, Wave Plus®	Woods Magnifier, Wave Plus®
Dental	Super Exam 50® LED, CoolSpot II®, Woods Magnifier	Super Exam 50® LED, CoolSpot II®, Woods Magnifier

Philips Burton offers Veterinary lights for nearly all applications. Utilize the table above to select the right light to meet your needs.

### AIM LED

- 45,000 lux (4181 fc) at 1 meter
- Saves 70% of the energy used by traditional sources, consuming only 56 watts of power
- 3 LED modules with a 20,000-hour life
- Adjustable focus using center handle



### AIM-50®

- 45,000 lux (4181 fc) at 1 meter
- Extremely smooth focus mechanism with center sterilizable handle
- 360° limitless arm-and-mounting-system-rotation around vertical axes (single ceiling)



### AIM-100®

- 63,000 lux (5852 fc) at 1 meter
- Uses three quartz halogen bulbs – each rated with a 2000-hour bulb life
- Compact 20" (51 cm) diameter lighthead
- Cool operation for long hours of use



### COOLSPOT II®

- Up to 60,000 lux (5574 fc) at 1 meter
- Bi-level dimmer for variable intensity
- Variable spot size of 2" - 10" (5 cm-25 cm)
- Focuses with central removable, autoclavable SteriHandle™



### OUTPATIENT® II

- 86,000 lux (7990 fc) at 24" (61 cm)
- 2000-hour average bulb life
- Focuses with central removable, autoclavable SteriHandle™
- Cool operation assured with heat filters and dichroic coatings



### SUPER EXAM 50® LED

- 39,000 lux (3623 fc) at 24" (61 cm)
- LED light module with 40,000-hour life
- Delivers value with longer life and lower energy consumption
- Drift-free K-arm ensures easy positioning



### LE-35

- 30,000 lux (2787 fc) at 18" (46 cm)
- Two to five times more bright white light than similar lights in its class
- 4000-hour average lamp life
- Drift-free K-arm with 42" (107 cm)



### GLEAMER®

- 45,000 lux (4181 fc) at 18" (46 cm)
- Cool housing temperature and cool light beam for hours of use without causing discomfort to physician or patient
- High-quality, space-saving 4-caster base



### AIM-200® OR

- 90,000 lux (8361 fc) at 1 meter
- Larger and adjustable light pattern
- Unlimited rotation and maximum reach of 63" (160 cm)
- Light field diameter: 10"-12" (25 cm-30 cm)



### APEX

- 120,000 lux (11,148 fc) at 1 meter
- Durable arm system recognized for drift-free, trouble-free operation
- Unlimited arm rotation
- Continuous intensity adjustment



### WOODS MAGNIFIER (COMING SOON)

- UV-A Ultraviolet light and white light with white light dimming
- LED light source provides long life
- Cordless operation (i.e. rechargeable battery powered)
- 5 diopter magnifying lens



### WAVE PLUS®

- Superior optics; 3.5 diopter fluorescent magnifier with 13" (33 cm) focal length
- Extremely bright; 2 independent 13-watt fluorescent bulb sources with 10,000-hour average bulb life
- 41" (105 cm), no-drift, articulating K-Arm



An industry-leading 5-year warranty\* means Philips Burton lights are sure to exceed your expectations!

\* Excludes Wave® Plus and LE-35

# Ordering Information

AIM-100®	
A100FL	AIM-100® Floorstand
A100SC	AIM-100® Single, Ceiling Mount
A100DC	AIM-100® Dual, Ceiling Mount
A100W	AIM-100® Wall Mount
AIM-200® OR	
A200FL	AIM-200® OR Floorstand
A200SC	AIM-200® OR Single, Ceiling Mount
A200DC	AIM-200® OR Dual, Ceiling Mount
A200W	AIM-200® OR Wall Mount
AIM-50®	
A50FL	AIM-50® Floorstand
A50SC	AIM-50® Single, Ceiling Mount
A50DC	AIM-50® Dual, Ceiling Mount
A50W	AIM-50® Wall Mount
AIM <sub>LED</sub>	
ALEDFL	AIM <sub>LED</sub> Floorstand
ALEDSC	AIM <sub>LED</sub> Single, Ceiling Mount
ALEDDC	AIM <sub>LED</sub> Dual, Ceiling Mount
ALEDW	AIM <sub>LED</sub> Wall Mount
APEX	
APXSC8	APEX Single, Ceiling Mount, 8 ft.
APXDC8	APEX Dual, Ceiling Mount, 8 ft.
APXSC10	APEX Single, Ceiling Mount, 10 ft.
APXDC10	APEX Dual, Ceiling Mount, 10 ft.
COOLSPOT II® (complete models only)	
CS316FL	CoolSpot II® Floorstand
CS316SC	CoolSpot II® Single, Ceiling Mount
CS316DC	CoolSpot II® Dual, Ceiling Mount
CS316W	CoolSpot II® Wall Mount
CS316ST	Single CoolSpot II® on Fastrac Single Head and Single Trolley
CS316DT	Dual CoolSpot II® on Fastrac Dual Head and Single Trolley

GLEAMER®	
GL30FL	Gleamer® Floorstand (with casters)
GL30W	Gleamer® with Extension Arm and Universal Wall Mount
GLW30FL	Gleamer® Wide Beam Spot with Floorstand (with casters)
GLW30W	Gleamer® Wide Beam Spot with Extension Arm and Universal Wall Mount
LE-35	
LE-35FL	LE-35 Floorstand
LE-35T	LE-35 Table Mount
OUTPATIENT® II (complete models only)	
OP216FL	Outpatient® II Floorstand
OP216SC	Outpatient® II Single, Ceiling Mount
OP216DC	Outpatient® II Dual, Ceiling Mount
OP216W	Outpatient® II Wall Mount
OP216ST	Single Outpatient® II on Fastrac Single Head and Single Trolley
OP216DT	Dual Outpatient® II on Fastrac Dual Head and Single Trolley
SUPER EXAM 50® <sub>LED</sub>	
SE50LFL	Super Exam 50® <sub>LED</sub> Floorstand
SE50LSC	Super Exam 50® <sub>LED</sub> Single, Ceiling Mount
SE50LW	Super Exam 50® <sub>LED</sub> Wall Mount
SE50LR	Super Exam 50® <sub>LED</sub> Rail Mount
SE50LT	Super Exam 50® <sub>LED</sub> Table Mount
WAVE® PLUS	
WM50FL	Wave® Plus Magnifier Floorstand
WM50W	Wave® Plus Magnifier with Wall/Table Mount
WOODS MAGNIFIER	
MAGWDSHH	Woods Magnifier, Hand Held
MAGWDSAM	Woods Magnifier, Arm Mount

For a complete listing of our products visit [burtonmedical.com](http://burtonmedical.com)

*NOTE: 120V standard; Models also available in 230V 50/60 Hz.  
Contact Philips Burton for additional information.*



Philips Burton  
21100 Lassen Street Chatsworth, CA 91311  
Phone: 818-701-8700 Fax: 818-701-8725  
[www.burtonmedical.com](http://www.burtonmedical.com)

© 2011 Philips Burton. All rights reserved. Product designs protected by copyright.  
We reserve the right to change details of design, materials and finishes.

LIT-0017, Rev. 5