



Major Surgery Lights

XenaLux™

Advanced design ensures top performance and a noticeably brighter light beam.

Powered by a 300-watt Xenon source and with light delivery via multi-fiber cables, Burton's XenaLux™ OR headlamp system offers superior intensity as well as remarkable reliability and value. The headlamp's optical system has a full range of articulation with a variable spot range, removable joystick adjustment and a new and improved vertically adjustable linkage arm for added reach and greater maneuverability. Unique design and optics allow the XenaLux™ Illuminator to provide superior features.

The Right Light

PHILIPS
burton

Headlamp Features

- Extremely bright light at working distance
- Variable spot with flawless circular beam
- Very light weight right angle headlamp
- Preferred for comfort
- Compact optical system
- Fully-adjustable, uniquely-designed linkage arm allows for vertical adjustment of the optical module
- Adapts to any Wolf port
- Attractive standard 3.5 bundle cables

Illuminator Features

- Produces highest intensity white light of any 300-watt Xenon source
- Quiet, cool venting system
- Rotary Turret for Wolf, Storz, Olympus, ACMI and Pilling connectors
- Hours-in-use meter
- Optional mobile floorstand offers vertical adjustment on upright assembly; the illuminator's shelf includes a handle for easy positioning




Summary of technical data

Illuminance:	120,000+ lux (11,148+ fc)
Bulb life:	500 hours
Power (light bulbs):	300-watt
Type of bulb:	Xenon
Product warranty:	3-year limited warranty
Origin:	Assembled in USA

Note: The optical data are nominal values based on measurements calculated according to UL/IEC performance standards at 120V



Sample Model No: XN10C = Complete XenaLux™ Headlight System



Accessories

#0880300

Mobile Floorstand for Illuminator



Philips Burton
21100 Lassen Street Chatsworth, CA 91311
Phone: 818-701-8700 Fax: 818-701-8725
www.burtonmedical.com

© 2011 Philips Burton. All rights reserved. Product designs protected by copyright.
We reserve the right to change details of design, materials and finishes.

LIT-0012, Rev. 2