



Dermatology/ Plastic Surgery Lighting

When precision counts, count on Burton.



burton®

Dermatology/ Plastic Surgery Lighting

Focusable light with shadow-free positioning provides precise illumination.

With extended, no-drift arms, Burton lights are reliable and easy to maneuver. High-intensity white light and true color rendering provide accurate visibility for dermatologists and plastic surgeons, while cool operation ensures patient comfort. UV lights and LED magnifiers with large viewing areas offer enhanced optical performance.



Burton proudly stands behind the reliability of its lighting products with an industry-leading 5-year limited warranty.*

Dept./Specialty	Good	Better	Best
Office			
Plastic	Outpatient II / Outpatient LED	AIM LED / ALED070	AIM-100 / ALED070
Derm	Super Bright Spot LED	Outpatient II / Outpatient LED	AIM LED / ALED070
Magnifier	UV Light	Wave LED	Epic LED Magnifier
ASC			
Plastic	AIM LED / AIM-50 / ALED070	AIM-100 / ALED070	ALED100
G OR	AIM-100 / ALED070	AIM-200 OR / ALED100	ALED100
Magnifier	UV Light	Wave LED	Epic LED Magnifier

Examination Lights

AIM-50®

- 45,000 lux (4181 fc) at 39.37" (1 m)
- 3100 K color temperature
- CRI of 96
- Three 50-watt quartz halogen bulbs
- 2,000-hour average rated bulb life



AIM® LED

- 45,000 lux (4181 fc) at 39.37" (1 m)
- 4300 K color temperature
- CRI of 92
- 50,000-hour average rated bulb life



AIM® HI LED

- AIM HI 100 - 100,000 lux (9290 fc) at 39.37" (1 m)
- AIM HI 70 - 70,000 lux (6503 fc) at 39.37" (1 m)
- 4300 K color temperature
- CRI of >95
- 50,000 hours average rated bulb life



OUTPATIENT® II

- 86,000 lux (7990 fc) at 24" (61 cm)
- 3300 K color temperature
- CRI of 94
- Three 50-watt halogen bulbs
- 2,000-hour average rated bulb life



OUTPATIENT® LED

- 90,000 lux (8361 fc) at 24" (61 cm)
- 3700 K color temperature
- CRI of 95
- High R9>87 to accurately render skin tones
- 50,000-hour average rated bulb life



SUPER BRIGHT SPOT® LED

- 20,000 lux (1858 fc) at 16" (40.64 cm)
- 5200 K color temperature
- CRI of >90
- 50,000-hour average rated bulb life



Surgical Lights

AIM-100®

- 63,000 lux (5852 fc) at 39.37" (1 m)
- 3500 K color temperature
- CRI of 97
- Three 35-watt quartz halogen bulbs
- 2,000-hour average rated bulb life



AIM-200® OR

- 90,000 lux (8161 fc) at 39.37" (1 m)
- 3700 K color temperature
- CRI of 96
- Larger and adjustable light pattern



APEX

- 120,000 lux (11,148 fc) at 39.37" (1 m)
- 4200 K color temperature
- CRI of 95
- Unlimited arm rotation



Magnifiers

UV LIGHT

- 5 diopter magnifier (2.25x magnification)
- Focal length of 8" (20 cm)
- Fluorescent bulbs produce UVA at approximately 365 nm wavelength
- 6,000-hour average rated bulb life



WAVE LED

- 4,600 lux (427 fc) at 11" (28 cm)
- 4000 K color temperature
- CRI of 80
- 50,000-hour average rated bulb life



EPIC® LED MAGNIFIER

- 3.5 diopter magnifier (1.88x magnification)
- 1,360 lux (126 fc) at 8" (20 cm)
- 3,830 lux (356 fc) at 11" (28 cm)
- 4000 K color temperature
- CRI of 90
- 5,000-hour average rated bulb life



Find additional light products
at burtonmedical.com

Ordering Information

AIM-50	
A50FL	AIM-50 Floor Stand
A50SC	AIM-50 Single Ceiling Mount
A50DC	AIM-50 Double Ceiling Mount
A50W	AIM-50 Wall Mount
AIM LED	
ALEDFL	AIM LED Floor Stand
ALEDSC	AIM LED Single Ceiling Mount
ALEDDC	AIM LED Double Ceiling Mount
ALEDW	AIM LED Wall Mount
AIM HI LED	
ALED070FL	ALED070 Floor Stand
ALED070SC	ALED070 Single Ceiling
ALED070DC	ALED070 Dual Ceiling
ALED070W	ALED070 Wall
ALED100FL	ALED100 Floor Stand
ALED100SC	ALED100 Single Ceiling
ALED100DC	ALED100 Dual Ceiling
ALED100W	ALED100 Wall
OUTPATIENT II (COMPLETE MODELS ONLY)	
OP216FL	Outpatient II Floor Stand
OP216SC	Outpatient II Single Ceiling Mount
OP216DC	Outpatient II Double Ceiling Mount
OP216W	Outpatient II Wall Mount
OP216ST	Outpatient II Single Fastrac Mount
OP216DT	Outpatient II Double Fastrac Mount
OUTPATIENT LED	
OPLEDFL	Outpatient LED Floor Stand
OPLEDSC	Outpatient LED Single Ceiling Mount
OPLEDDC	Outpatient LED Double Ceiling Mount
OPLEDW	Outpatient LED Wall Mount
SUPER BRIGHT SPOT LED	
SBL-FL	Super Bright Spot® LED w/ Floor Stand (120V Wall Plug)
SBL-T	Super Bright Spot® LED w/ Table Mount (120V Wall Plug)
SBL-W	Super Bright Spot® LED w/ Wall Mount (120V Wall Plug)
SBL-P	Super Bright Spot® LED w/ Pole Mount (120V Wall Plug)
SBL-FL-INT	Super Bright Spot® LED w/ Floor Stand with International Wall Plug Kit
SBL-T-INT	Super Bright Spot® LED w/ Table Mount with International Wall Plug Kit
SBL-W-INT	Super Bright Spot® LED w/ Wall Mount with International Wall Plug Kit
SBL-P-INT	Super Bright Spot® LED w/ Pole Mount with International Wall Plug Kit

NOTE: 120V standard; Additional voltage models also available. Contact Burton for additional information.



AIM-100	
A100FL	AIM-100 Floor Stand
A100SC	AIM-100 Single Ceiling Mount
A100DC	AIM-100 Double Ceiling Mount
A100W	AIM-100 Wall Mount
AIM-200 OR	
A200FL	AIM-200 Floor Stand
A200SC	AIM-200 Single Ceiling Mount
A200DC	AIM-200 Double Ceiling Mount
A200W	AIM-200 Wall Mount
APEX	
APXSC8	APEX 8' Single Ceiling Mount
APXSC10	APEX 10' Single Ceiling Mount
APXDC8	APEX 8' Double Ceiling Mount
APXDC10	APEX 10' Double Ceiling Mount
UV LIGHT	
UV501	UV Light, 2 UV Bulbs without Magnifier
UV502	UV Light, 2 UV Bulbs with Magnifier
UV503	UV Light, 2 UV, 2 White Bulbs with Magnifier
WAVE LED	
WML35W45EC120	Wave LED (standard table edge clamp included)
Optional Accessory:	
WMLWFS	Weighted Floor Stand
EPIC LED MAGNIFIER	
EMLFL	Epic LED Magnifier, Floor Stand, NEMA 5-15P
EMLT	Epic LED Magnifier, Table Mount, NEMA 5-15P
EMLW	Epic LED Magnifier, Wall Mount, NEMA 5-15P
EMLFL01	Epic LED Magnifier, Floor Mount, CEE 7/7
EMLFL02	Epic LED Magnifier, Floor Mount, AS/NZS 3112
EMLFL03	Epic LED Magnifier, Floor Mount, BSI363
EMLT01	Epic LED Magnifier, Table Mount, CEE 7/7
EMLT02	Epic LED Magnifier, Table Mount, AS/NZS 3112
EMLT03	Epic LED Magnifier, Table Mount, BSI363
EMLW01	Epic LED Magnifier, Wall Mount, CEE 7/7
EMLW02	Epic LED Magnifier, Wall Mount, AS/NZS 3112
EMLW03	Epic LED Magnifier, Wall Mount, BSI363

* To register your Burton light(s), please visit:
<http://www.burtonmedical.com/warranty-registration/>

Find additional light products and accessories at burtonmedical.com

© 2018 Burton Medical. All rights reserved.
Burton Medical reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

burton®

Burton Medical, LLC
2300 W Windsor Ct., Suite C
Addison, IL 60101 USA
Phone (800) 444-9909, Fax (800) 765-1770
Burtonmedical.com