



OB/GYN Lighting

Outstanding performance for the practitioner,
exceptional comfort for the patient.



burton®

OB/GYN Lighting

High-intensity illumination with minimal heat output creates the perfect light for both the physician and patient.

Burton lights deliver precisely controlled natural white light that is so important for an accurate gynecological exam. Whether you're looking for a pre-focused spot or a flood-type light, Burton offers a variety of lights to meet your need. Compact sizing and no-drift arms make these lights easy to work with and perfect for any facility. Our LED models provide overall cost of ownership advantages with low maintenance and energy savings.



Burton proudly stands behind the reliability of its lighting products with an industry-leading 5-year limited warranty.*

Dept./Specialty	Good	Better	Best
Office			
OB/GYN	Gleamer / Super Bright Spot	Nova Exam LED / Super Bright Spot LED	CoolSpot II
Hospital			
LDR	Outpatient II / Outpatient LED	CoolSpot II	AIM LED / ALED070 / ALED100
Exam	Gleamer, Super Bright Spot	Nova Exam LED	Outpatient II / Outpatient LED Exam Light / ALED70 / ALED100

Examination Lights

AIM® HI LED

- AIM HI 100 - 100,000 lux (9290 fc) at 39.37" (1 m)
- AIM HI 70 - 70,000 lux (6503 fc) at 39.37" (1 m)
- 4300 K color temperature
- CRI of >95
- 50,000 hours average rated bulb life



GLEAMER®

- 45,000 lux (4181 fc) at 18" (46 cm)
- 3200 K color temperature
- CRI of 91
- 35-watt halogen bulb
- 2,000-hour average rated bulb life



NOVA EXAM® LED

- 25,000 lux (low level) (2323 fc) at 18" (46 cm)
- 50,000 lux (high level) (4645 fc) at 18" (46 cm)
- 4600 K color temperature
- CRI of 93
- 40,000-hour average rated bulb life



COOLSPOT® II

- Up to 60,000 lux (5574 fc) at 39.37" (1 m)
- 4300 K color temperature
- CRI of 94
- 600-hour average rated bulb life



OUTPATIENT® II

- 86,000 lux (7990 fc) at 24" (61 cm)
- 3300 K color temperature
- CRI of 94
- Three 50-watt halogen bulbs
- 2,000-hour average rated bulb life



OUTPATIENT LED®

- 90,000 lux (8361 fc) at 24" (61 cm)
- 3700 K color temperature
- CRI of 95
- High R9>87 to accurately render skin tones
- 20,000-hour average rated bulb life



OUTPATIENT® II/COOLSPOT® II COMBO

- Outpatient II:
 - 86,000 lux (7990 fc) at 24" (61 cm)
 - 3300 K color temperature
- CoolSpot II:
 - Up to 60,000 lux (5574 fc) at 39.37" (1 m)
 - 4300 K color temperature



AIM LED

- 45,000 lux (4181 fc) at 39.37" (1 m)
- 4300 K color temperature
- CRI of 92
- 50,000-hour average rated bulb life



SUPER BRIGHT SPOT® LED

- 20,000 lux (1858 fc) at 16" (40.64 cm)
- 5200 K color temperature
- CRI of >90
- 50,000-hour average rated bulb life



SUPER BRIGHT SPOT

- 56,000 lux (5185 fc) at 18" (46 cm)
- 3400 K color temperature
- CRI of 98
- 20-watt quartz halogen bulb
- 2,000-hour average rated bulb life



Find additional light products
at burtonmedical.com

Ordering Information

AIM HI LED	
ALED070FL	ALED070 Floor Stand
ALED070SC	ALED070 Single Ceiling
ALED070DC	ALED070 Dual Ceiling
ALED070W	ALED070 Wall
ALED100FL	ALED100 Floor Stand
ALED100SC	ALED100 Single Ceiling
ALED100DC	ALED100 Dual Ceiling
ALED100W	ALED100 Wall

GLEAMER	
GL30FL	Gleamer Floor Stand
GL30W	Gleamer Wall Mount

NOVA EXAM LED	
NXFL	Nova Exam LED Floor Stand
NXT	Nova Exam LED Table mount
NXW	Nova Exam LED Wall Mount
NXC	Nova Exam LED Ceiling Mount

COOLSPOT II (complete models only)	
CS316FL	CoolSpot II Floor Stand
CS316SC	CoolSpot II Single Ceiling Mount
CS316DC	CoolSpot II Double Ceiling Mount
CS316W	CoolSpot II Wall Mount
CS316ST	CoolSpot II Single Fastrac mount
CS316DT	CoolSpot II Double Fastrac mount

OUTPATIENT II (complete models only)	
OP216FL	Outpatient II Floor Stand
OP216SC	Outpatient II Single Ceiling Mount
OP216DC	Outpatient II Double Ceiling Mount
OP216W	Outpatient II Wall Mount
OP216ST	Outpatient II Single Fastrac Mount
OP216DT	Outpatient II Double Fastrac Mount

OUTPATIENT LED	
OPLEDFL	Outpatient LED Exam Light Floor Stand
OPLEDSC	Outpatient LED Exam Light Single Ceiling Mount
OPLEDDC	Outpatient LED Exam Light Double Ceiling Mount
OPLEDW	Outpatient LED Exam Light Wall Mount
OPLEDST	Outpatient LED Exam Light Single Fastrac Mount
OPLEDDT	Outpatient LED Exam Light Double Fastrac Mount

OUTPATIENT II / COOLSPOT II COMBO	
OC516CC	Outpatient II / CoolSpot II Ceiling Mount
OC516CT	Outpatient II / CoolSpot II Fastrac Mount

AIM LED	
ALEDFL	AIM LED Floor Stand
ALEDSC	AIM LED Single Ceiling Mount
ALEDDC	AIM LED Double Ceiling Mount
ALEDW	AIM LED Wall Mount

SUPER BRIGHT SPOT LED	
SBL-FL	Super Bright Spot® LED w/ Floor Stand (120V Wall Plug)
SBL-T	Super Bright Spot® LED w/ Table Mount (120V Wall Plug)
SBL-W	Super Bright Spot® LED w/ Wall Mount (120V Wall Plug)
SBL-P	Super Bright Spot® LED w/ Pole Mount (120V Wall Plug)
SBL-FL-INT	Super Bright Spot® LED w/ Floor Stand with International Wall Plug Kit
SBL-T-INT	Super Bright Spot® LED w/ Table Mount with International Wall Plug Kit
SBL-W-INT	Super Bright Spot® LED w/ Wall Mount with International Wall Plug Kit
SBL-P-INT	Super Bright Spot® LED w/ Pole Mount with International Wall Plug Kit

SUPER BRIGHT SPOT	
SB40FL	Super Bright Spot Floor Stand
SB40W	Super Bright Spot Wall Mount

NOTE: 120V standard; Additional voltage models also available. Contact Burton for additional information.

* To register your Burton light(s), please visit: <http://www.burtonmedical.com/warranty-registration/>

Find additional light products and accessories at burtonmedical.com



© 2018 Burton Medical. All rights reserved. Burton Medical reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

burton®

Burton Medical, LLC
2300 W Windsor Ct., Suite C
Addison, IL 60101 USA
Phone (800) 444-9909, Fax (800) 765-1770
Burtonmedical.com