



Task Light

LE 50[®] LED

Burton's economical multi-purpose LED light.

Burton's LE 50 LED is a cost-effective, flex-arm, medical task, and reading light. Its illumination far exceeds the competition.



burton[®]

Applications

Medical Task

- Ambulatory Surgery Centers
- Urgent Care
- Emergency / Triage Centers
- Long Term Care
- Veterinary
- School Nurse


Patient Reading

- Patient Rooms
- Senior Living / Long Term Care

Features

- High-quality illumination
- Durable construction
- Energy efficient, low maintenance LED
- Easy positioning light head
- Convenient on/off switch
- Standard or wide beam models available
- Mobile floor stand, or space-saving wall or table mount models available

LE 50® LED product dimensions



Summary of technical data

LE 50® LED

Illuminance:	Standard Beam: 40,000 Lux at 18"; 75,000 Lux at 12" Wide Beam: 9,000 Lux at 18"; 19,000 Lux at 12"
Color Temperature:	4100 Kelvin
CRI (Color Rendering Index):	95
Power Supply:	Universal Voltage 100/230V External Adapter
Power:	7 Watts
Light Source:	Light Emitting Diode (LED)
Rated Life of LED:	50,000-hour average rated LED life
Mounting Options:	Floor Stand, Wall, Table, Rail
Arm:	Flexible Gooseneck
Certifications/Approvals:	UL / cUL
Product Warranty:	2-year limited warranty
Origin:	Assembled in USA

Ordering Information

LE 50® LED - STANDARD BEAM



Universal Input 100V/230V	Description
LE50FL	LE 50 LED Standard Beam Mobile Floor Stand
LE50W	LE 50 LED Standard Beam Extension Arm & Universal Wall Mount
LE50T	LE 50 LED Standard Beam Table Mount

LE 50® LED - WIDE BEAM

LEW50FL	LE 50 LED Wide Beam Mobile Floor Stand
LEW50W	LE 50 LED Wide Beam Extension Arm & Universal Wall Mount
LEW50T	LE 50 LED Wide Beam Table Mount

ACCESSORIES

Model	Description
4000068	Din Style Rail Mount 25mm X 10mm
4000067	Fairfield Style Rail Mount
6000431	Wall Mount Extension 15.5"



burton®

Burton Medical, LLC
2300 W Windsor Ct., Suite C
Addison, IL 60101 USA
Phone (800) 444-9909, Fax (800) 765-1770
Burtonmedical.com

© 2021 Burton Medical. All rights reserved. Burton Medical reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.