



Task Light

Gleamer LED[®]

Burton's all-purpose LED light.

If your facility is searching for an all-purpose, rugged, compact LED light, Burton's Gleamer LED[®] is your solution. The Gleamer LED[®] designed to meet a healthcare facilities' most demanding specifications for light-intensity, cool operation, compact design, and value. With an impressive 50,000 lux output, a sturdy gooseneck arm for precise light head positioning and stability. Built with superior-quality components, the Gleamer LED[®] ensures exceptional reliability – long life, maximum functionality, and outstanding brightness for a wide range of applications.

Features

- Cool housing temperature and cool light beam; light can be in continuous use for hours without causing discomfort to practitioner or patient
- Convenient positioning with ergonomic positioning ridges on light head
- High-quality, space-efficient 5-caster base designed for smooth and flexible mobility. Casters are lockable for maximum stability and secure positioning.



burton[®]

Summary of technical data

GLEAMER LED

Illuminance:	Narrow Beam: 50,000 Lux at 18" / 460 mm Wide Beam: 9,000 Lux at 18" / 460 mm
Spot Size:	Narrow Beam: 165 mm (6.5 inch) at 18" / 460 mm Wide Beam: 420 mm (16.4 inch) at 18" / 460 mm
Color Temperature:	3,500 Kelvin
Light Source:	LED
Rated Life:	50,000-hour average bulb life
Power:	7 Watts
Power Supply:	100-240V AC 50/60Hz
Certifications / Approvals:	UL153/ CSA250.4
Product Warranty:	2-year limited warranty

Ordering Information


GLEAMER LED

Narrow Beam	Description
GLEDFL	Narrow Beam Floor
GLEDW	Narrow Beam Wall / Exam Table with Extension Arm
GLEDT	Narrow Beam Table
Wide Beam	Description
GWLEDFL	Wide Beam Floor
GWLEDW	Wide Beam Wall / Exam Table with Extension Arm
GWLEDT	Wide Beam Table

Gleamer LED® product dimensions



Mounting Holes fit for Exam Table / Procedure Chair Brackets



EnabLED
Licensing Program for LED
Luminaires and Retrofit Bulbs

Member
Burton Medical, LLC



© 2024 Burton Medical. All rights reserved. Burton Medical reserves the right to make changes in specifications and/or to discontinue any product at any time without notice

burton®

Burton Medical, LLC
2300 W. Windstor Ct., Suite C
Addison, IL 60101 USA
Phone (800) 444-9909, Fax (800) 765-1770