



Examination Light

Outpatient II / CoolSpot II Combination Models

The minimal heat output and high color temperature of the Outpatient II along with the incredibly pure, white light of the CoolSpot II combined. Available in Ceiling Mount or Fastrac Mount.

Outpatient II Features:


- 86,000 Lux at 24" (61 cm)
- 3300 K Color Temperature
- CRI (Color Rendering Index) of 94
- Focuses with central removable, autoclavable SteriHandle

CoolSpot II Features:


- 60,000 Lux at 1 meter
- 4300 K Color Temperature
- CRI (Color Rendering Index) of 94
- Bi-level dimmer for variable intensity



burton®



OC516CC



OC516CT

Summary of technical data

OUTPATIENT II

Illuminance:	86,000 Lux at 24" (61 cm)
Color Temperature:	3300 Kelvin
CRI (Color Rendering Index):	94
Light Source:	Halogen
Rated Life:	2000-hour average bulb life
Power:	50 Watts
Power:	Ceiling Mount or Fastrac Mount
Product Warranty:	5 Year Limited Warranty

Summary of technical data

COOLSPOT II

Illuminance:	60,000 Lux at 1 meter
Color Temperature:	4300 Kelvin
CRI (Color Rendering Index):	94
Light Source:	Halogen
Rated Life:	600-hour average bulb life
Power:	150 Watts
Power:	Ceiling Mount or Fastrac Mount
Product Warranty:	5 Year Limited Warranty

Ordering Information

OUTPATIENT II / COOLSPOT II COMBINATION MODEL

120v Model	230/240v Model	230/240v SP Model*	100v Model	Description
OC516CC	OC525CC	O202CC25	OC515CC	Combo Ceiling Mount
OC516CT	OC525CT		OC515CT	Combo Trolley Mount

*SP Model refers to Spanish Language Manuals.

ACCESSORIES

Model No.	Description
0006130PK	Outpatient II Replacement Bulbs, 3/Box
1008465	Outpatient II Handle Replacement
0007006PK	CoolSpot II Replacement Bulbs, 3/Box
1008466B	CoolSpot II Handle Replacement
0008100PK	Light Handle Cover, Sterile, 25/Box
1008259	42" Drop Tube
1008060	Special Length, Max. 58" Drop Tube



© 2022 Burton Medical. All rights reserved. Burton Medical reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

burton®

Burton Medical, LLC
 2300 W Windsor Ct., Suite C
 Addison, IL 60101 USA
 Phone (800) 444-9909, Fax (800) 765-1770
 Burtonmedical.com