



Task Light

LE 50[®] LED

Burton's economical multi-purpose LED light.

Burton's LE 50 LED is a cost-effective, flex-arm, medical task, and reading light. Its illumination far exceeds the competition.



burton[®]

Applications

Medical Task

- Ambulatory Surgery Centers
- Urgent Care
- Emergency / Triage Centers
- Long Term Care
- Veterinary
- School Nurse

Patient Reading

- Patient Rooms
- Senior Living / Long Term Care

Features


- High-quality illumination
- Durable construction
- Energy efficient, low maintenance LED
- Easy positioning light head
- Convenient on/off switch
- Mobile floor stand, or space-saving wall or table mount models available

Summary of technical data

LE 50® LED

Illuminance:	Standard Beam: 40,000 Lux at 18"; 75,000 Lux at 12"
Color Temperature:	4100 Kelvin
CRI (Color Rendering Index):	95
Power Supply:	Universal Voltage 100/230V External Adapter
Power:	7 Watts
Light Source:	Light Emitting Diode (LED)
Rated Life of LED:	50,000-hour average rated LED life
Mounting Options:	Floor Stand, Wall, Table
Arm:	Flexible Gooseneck
Certifications/Approvals:	UL / cUL
Product Warranty:	1-year limited warranty
Origin:	Assembled in USA

LE 50® LED product dimensions



Ordering Information

LE 50® LED - STANDARD BEAM

Universal Input 100V/230V	Description
LE50FL	LE 50 LED Standard Beam Mobile Floor Stand
LE50W	LE 50 LED Standard Beam Universal Wall Mount
LE50T	LE 50 LED Standard Beam Table Mount

ACCESSORIES

Model	Description
5000167	Wall Bracket

5000167



Mounting holes on 4 corners

burton®

Burton Medical, LLC
2300 W Windsor Ct, Suite C
Addison, IL 60101 USA
Phone (800) 444-9909, Fax (800) 765-1770
Burtonmedical.com

© 2021 Burton Medical. All rights reserved. Burton Medical reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.