

Headlights

LED HL70

Comfortable, lightweight, high quality headlight that helps you go the distance.

Features:

- 70,000 lux at 16" (40.64 cm) working distance
- LED Head Light module provides exceptional light quality
- Long Battery Life 2 hours 45 minutes at full light intensity, 3.5 hour recharge
- Variable spot size .78" - 4.33"
- Variable intensity control 0 - 70,000 lux at 16" (40 cm)
- Lightweight Head Lamp - 9 oz (225 g)



LED HL70



burton®

Summary of technical data

LED HEADLIGHT HL70

Illuminance:	70,000 Lux at 16" (40.64 cm)
Color Temperature:	4,500 Kelvin
CRI (Color Rendering Index):	75
Light Source:	LED (Light Emitting Diode)
Rated Life of LED:	50,000 hours
Power:	10 Watts
Headband and LED Module Weight:	9 ounces (225 grams)
Certifications/Approvals	FDA Class I
Product Warranty:	3-year Limited Headband and LED Module 1-year Limited Battery

Ordering Information

LED HEADLIGHT HL70

Model	Description
HL70	Complete System - LED Headlight, Carrying Case, 2 Batteries, Battery Charger, 2 Joysticks Non-Sterile

ACCESSORIES

Model	Description
4000445	Battery Charger for HL70, 120-240V
4000442	Encrypted Battery Replacement HL70 Headlight
4000446	Joysticks Non-Sterile Autoclavable, Pack of 5



Applications:

- Surgery
- ENT
- Dermatology
- Plastic Surgery
- Veterinary
- Dental

